

TM

# DILATOR BALLOONS

ESOPHAGEAL, PYLORIC/ COLONIC, BILLIARY, ACHALASIA



# DILATOR BALLOONS™

ESOPHAGEAL, PYLORIC /COLONIC,BILLIARY, ACHALASIA

## Product Features : -

- Transparent Material
- Metal Marks For Easy Positioning
- Guide Wire Compatible
- Non Compliant Balloon
- Fast Inflation & Deflation

Catalog No.	Type	Balloon OD	Balloon Length	Catheter OD	Catheter Length
EBD - 0110080	Esophageal	10 mm	8 cm	2.4 mm	180 cm
EBD - 0112080	Esophageal	12 mm	8 cm	2.4 mm	180 cm
EBD - 0114080	Esophageal	14 mm	8 cm	2.4 mm	180 cm
EBD - 0116080	Esophageal	16 mm	8 cm	2.4 mm	180 cm
EBD - 0118080	Esophageal	18 mm	8 cm	2.4mm	180 cm
EBD - 0210060	Pyloric / Colonic	10 mm	6 cm	2.4 mm	180 cm
EBD - 0212060	Pyloric /Colonic	12 mm	6 cm	2.4 mm	180 cm
EBD - 0214060	Pyloric / Colonic	14 mm	6 cm	2.4 mm	180 cm
EBD - 0216060	Pyloric / Colonic	16 mm	6 cm	2.4 mm	180 cm
EBD - 0218060	Pyloric / Colonic	18 mm	6 cm	2.4 mm	180 cm
EBD - 0220060	Pyloric / Colonic	20 mm	6 cm	2.4 mm	180 cm
EBD - 0306040	Biliary	6 mm	4 cm	2.4 mm	200 cm
EBD - 0308040	Biliary	8 mm	4 cm	2.4 mm	200 cm
EBD - 0310040	Biliary	10 mm	4 cm	2.4 mm	200 cm
EBD - 0430080	Achalasia	30 mm	8 cm	5 mm	90 cm
EBD - 0435080	Achalasia	35 mm	8 cm	5 mm	90 cm
EBD - 0440080	Achalasia	40 mm	8 cm	5 mm	90 cm

Manufactured by:

 **Endo-Med**<sup>®</sup>  
Technologies Pvt.Ltd.

New Delhi - INDIA  
Ph. :+ 91 11 22048222  
Fax :+ 91 11 22447748  
URL : www.endomed.in  
E-mail : info@endomed.in

Distributed by: