

# NITINOL GUIDEWIRE <sup>TM</sup>



# NITINOL GUIDEWIRE<sup>TM</sup>

## Product Features : -

- Nitinol Core.
- Highly Radiopaque Distal 5 c.m Tip
- PTFE Jacket
- Hydrophilic Coated
- Twin Color & Highly Lubricious

Catalog No.	Diameter	Length
EGW - 0125260	0.025"	260 cm
EGW - 0125450	0.025"	450 cm
EGW - 0135260	0.035"	260 cm
EGW - 0135450	0.035"	450 cm

Other size & configurations are available on request

Manufactured by:

 **Endo-Med<sup>®</sup>**  
**Technologies Pvt.Ltd.**

New Delhi - INDIA  
Ph. :+91 11 22048222  
Fax :+91 11 22 44 7748  
URL :www.endomed.in  
E-mail :info@endomed.in

Distributed by: