


Track-WireTM

Track-wire may be left in place during sphincterotomy.



Specially designed tip for tight Strictures and torturous anatomy.

Track-Wire™

Product Features : -

- Nitinol Core.
- Radiopaque Distal - Tungsten Filled Tip (5 c.m.).
- Highly Lubricious PTFE Jacket and Hydrophilic Coated.

Catalog No.	Diameter	Length
EGW - 0225260	0.025"	260 cm
EGW - 0225450	0.025"	450 cm
EGW - 0235260	0.035"	260 cm
EGW - 0235450	0.035"	450 cm

Manufactured by:

 **Endo-Med**[®]
Technologies Pvt.Ltd.

New Delhi - INDIA
Ph. : +91-11-42448449
Fax : +91-11-22447748
URL : www.endomed.in
E-mail : info@endomed.in

Distributed by: